

*Specification for Powder free latex Examination Gloves OP-PERFECT EXCELLENT*

<i>Color</i>	<i>white</i>												
<i>Shape</i>	<i>Straight Fingers</i>												
<i>Cuff</i>	<i>Beaded</i>												
<i>Powder free latex Examination Gloves</i>													
<i>PHYSICAL DIMENSIONS</i>													
<i>Size</i>	<i>Length (mm)</i>		<i>Palm Width (mm)</i>				<i>Thickness(mm)</i>						
	<i>Standard</i>	<i>Test</i>	<i>Standard</i>	<i>Test</i>	<i>Standard</i>								
					<i>Cuff</i>	<i>Palm</i>	<i>Fin</i>						
					<i>Smooth 0.08 Min</i>	<i>Thextured 0.11Min</i>	<i>Thextured 0.11Min</i>						
					<i>Test</i>								
<i>Small</i>	$\geq 230$	258-270	$80\pm 10$	84-87		0.08-0.09		0.11-0.12		0.12-0.13			
<i>Medium</i>		258-270	$95\pm 10$	96-99		0.08-0.09		0.11-0.12		0.12-0.13			
<i>Large</i>		258-270	$110\pm 10$	106-109		0.08-0.09		0.11-1.12		0.12-0.13			
<i>Tensile Properties</i>		<i>Standard</i>		<i>Test</i>									
		<i>EN</i>	<i>ISO</i>	<i>1</i>	<i>2</i>	<i>3</i>	<i>4</i>	<i>5</i>	<i>6</i>	<i>7</i>	<i>8</i>	<i>9</i>	<i>10</i>
<i>Before aging</i>	<i>Max force (N)</i>	3.6	7	7.8	8	7.9	7.7	7.8	7.6	7.9	7.7	7.8	7.9
	<i>Elongation (%)</i>	500	650	710	710	720	730	770	730	740	740	750	740
<i>After aging</i>	<i>Max force (N)</i>	3.6	6	6.9	6.8	7.1	7.2	6.8	6.9	7	7.1	6.9	6.8
	<i>Elongation (%)</i>	400	500	650	680	660	680	670	670	670	680	660	670
<i>Protein (<math>\mu\text{gr}/\text{gr}</math>)</i>		50	50	37.5	40	38.5	40.3	37.1	38.4	37.2	39.1	38.6	38.5
<i>Powder (mg)</i>		0.67	0.67	0.53	0.55	0.58	0.58	0.53	0.57	0.52	0.59	0.54	0.59