

## Specification for powder Free LATEX HARIR surgical gloves

<i>Color</i>	<i>Natural color</i>												
<i>Shape</i>	<i>Curve Fingers</i>												
<i>Cuff</i>	<i>Beaded</i>												
<i>Powder Free latex HARIR surgical Gloves</i>													
<b>PHYSICAL DIMENSIONS</b>													
<i>Size</i>	<i>Length (mm)</i>		<i>Palm Width (mm)</i>				<i>Thickness(mm)</i>						
	<i>Standard</i>	<i>Test</i>	<i>Standard</i>	<i>Test</i>	<i>Standard</i>								
					<i>Cuff</i>			<i>Palm</i>			<i>Fin</i>		
					<i>Smooth 0.10 Min</i>			<i>Thextured 0.13 Min</i>			<i>Thextured 0.13 Min</i>		
					<i>Test</i>								
6.5	≥260	285	83±5	85	0.14	0.15	0.22						
7	≥270	280	89±5	91	0.14	0.16	0.22						
7.5	≥270	284	95±5	99	0.14	0.15	0.21						
8	≥270	286	102±6	105	0.14	0.16	0.22						
<i>Tensile Properties</i>		<i>Standard</i>			<i>Test</i>								
		<i>EN</i>	<i>ISO</i>	<i>1</i>	<i>2</i>	<i>3</i>	<i>4</i>	<i>5</i>	<i>6</i>	<i>7</i>	<i>8</i>	<i>9</i>	<i>10</i>
<i>Before aging</i>	<i>Max force (N)</i>	9	12.5	17.2	17.1	16.9	17	17.2	16.8	16.9	17.2	17.1	17.2
	<i>Elongation (%)</i>	-	700	996	1010	1008	997	1006	1000	998	995	997	998
	<i>Initial Modulus (MPa)</i>	-	2	0.81	0.75	0.85	0.90	0.88	0.84	0.87	0.85	0.89	0.84
<i>After aging</i>	<i>Max force (N)</i>	9	9.5	15.5	15.8	16	15.9	15.8	16.1	15.8	15.7	16	15.8
	<i>Elongation (%)</i>	-	550	880	864	871	897	859	899	864	901	904	885
<i>Protein (µgr/gr)</i>		50	50	48	45	49	47	46	39	47	45	47	48
<i>Powder (mg)</i>		>2	0.6	0.58	0.55	0.57	0.55	0.54	0.50	0.56	0.49	0.58	0.53